




# PMCS-150

## Cooled High Speed PMT Module for TCSPC

- Fast TCSPC Instrument Response: < 120 ps FWHM**
- Internal Cooler: Low Dark Count Rate**
- Internal GHz Preamplifier: High Output Amplitude**
- Internal Overload Protection**
- Internal High Voltage Generation**
- Simple +12V Power Supply**
- Direct Interfacing to all bh Photon Counting Devices**
- Standard C Mount Adapter**
- Small Size: 40 x 40 x 120 mm**



The PMCS-150 is a cooled detector module for TCSPC applications. It contains a fast miniature PMT along with a Peltier cooler, a high voltage generator, a GHz pulse amplifier and a current sensing circuit. Due to the high gain and bandwidth of the device a single photon yields an output pulse with an amplitude in the range of 100 to 200 mV and a pulse width of 1.5 ns. The TCSPC instrument response function (IRF) has a width of less than 120 ps FWHM. It decreases to less than 100 ps for small sizes of the illuminated spot on the photocathode. The detector is operated with a fixed gain and from a single +12 V power supply. Overload conditions are internally detected by sensing the PMT output current. When overload occurs the PMT operating voltage is regulated down so that the output current remains on a safe level. Alternatively, the PMCS-150 can be operated from a bh DCU series detector controller module. In that case, the controller provides for power supply, gain control, and overload shutdown. The PMCS-150 interfaces directly to all bh SPC series time-correlated single photon (TCSPC) devices.







Becker & Hickl GmbH  
Nunsdorfer Ring 7-9  
12277 Berlin  
Tel. +49 / 30 / 212 800 20  
FAX +49 / 30 / 212 800 213  
<http://www.becker-hickl.com>  
email: [info@becker-hickl.com](mailto:info@becker-hickl.com)



# PMCS-150

	PMCS-150-04	PMCS-150-110	PMCS-150-113	PMCS-150-210	PMCS-150-20
Cathode type	Multialkali, UV	Super Bialkali	Super Bialkali, UV	Ultra Bialkali	Extended Red
Wavelength Range (nm)	185 to 870	230-700	185 to 700	230 to 700	300 to 900
Dark Counts (with cooling enabled <sup>1)</sup> )	40	20	20	40	500
Cathode Diameter			8 mm		
IRF width, 1mm spot / full cathode area, typ., FWHM			100 ps / 120 ps		
IRF shift with x-y position on cathode, see diagrams			50 ps		
IRF shift with magnetic field (150 Gauss), see diagram			10 ps		
Single-Electron Response Width			1.5 ns, FWHM, typ. value		
Single-Electron Response Amplitude			-100 to -200 mV		
Count Rate (Continuous)			> 5 MHz		
Count Rate (Peak, < 1 us)			> 100 MHz		
Internal Overload Protection			Down-regulation of PMT voltage		
Recovery from overload			Automatically, when light intensity has returned to normal		
Overload Indicator			LED		
Overload Signal			TTL / CMOS, open drain, active low		
Detector Signal Output Connector			SMA		
Output Impedance			50 Ω		
Power Supply			+ 12 V, 300 mA <sup>2)</sup>		
Dimensions (width x height x depth)			40 mm x 40 mm x 120 mm		
Optical Adapter			C-Mount female		
Fibre Coupling			SMA 905 or FC, on request		

- 1) Detector mounted on heat sink, temperature of detector base <25 °C  
 2) 100 mA when cooler is disabled



### 15 pin power-supply connector

- |         |                                       |        |                                |
|---------|---------------------------------------|--------|--------------------------------|
| 1, 2    | +12V                                  | 10, 11 | not used                       |
| 3       | /Disable Cooler, leave open to enable | 12     | not used                       |
| 4, 5    | GND                                   | 13     | Vgain, optional 0..1V from DCC |
| 6       | /Shutdown, leave open to enable       | 14     | not used                       |
| 7, 8, 9 | GND                                   | 15     | /OVLD, optional, to DCC        |

### International Sales Representatives



US:  
**Boston Electronics Corp**  
 tcspc@boselec.com  
 www.boselec.com



UK:  
**Photonic Solutions PLC**  
 sales@psplc.com  
 www.psplc.com



Japan:  
**Tokyo Instruments Inc.**  
 sales@tokyoinst.co.jp  
 www.tokyoinst.co.jp



China:  
**DynaSense Photonics Co. Ltd.**  
 info@dyna-sense.com  
 www.dyna-sense.com