

SOPHIA[®]-XO 2048B

SOPHIA[®]-XO 4096B

LARGE-FORMAT, DIRECT-DETECTION CCD CAMERAS FOR SOFT X-RAYS

Datasheet



TELEDYNE
PRINCETON INSTRUMENTS
Everywhere you look™

Breakthrough Soft X-ray Performance




The SOPHIA-XO 2048B and 4096B represent the latest generation of EUV/VUV and soft x-ray direct-detection cameras:

- ▶ Back-illuminated CCD sensors with >95% QE (~5 eV to 30 keV range)
- ▶ 2k x 2k and 4k x 4k formats; 13.5 and 15 micron pixels
- ▶ High frame rates with up to 4-port readout
- ▶ Cooling down to -90°C using liquid or air

Applications include:

VUV/EUV/XUV Imaging | X-ray Diffraction | X-ray Microscopy | X-ray Holography | X-ray Spectroscopy | X-ray Plasma


Princeton Instruments SOPHIA-XO soft x-ray cameras deliver speed, resolution, and sensitivity.



Designed for low-flux applications...

SOPHIA-XO delivers optimum performance:

- ▶ Low read noise
- ▶ High QE (>95% peak)
- ▶ Wide dynamic range with up to 18-bit readout



When speed is paramount...


SOPHIA-XO features the newest readout electronics:

- ▶ 1-, 2-, or 4-port simultaneous readout
- ▶ Multiple ADC speeds (up to 16 MHz)
- ▶ Binning and ROI readout
- ▶ Custom readout modes for microsecond exposures

We give you the most options!


Most complex experiments require flexibility... SOPHIA-XO has what you need:

- ▶ Air or liquid cooling
- ▶ LightField with built-in math engine
- ▶ Microsoft® Windows® 10 or Linux® 64-bit operating system support
- ▶ Compatibility with Python®, Linux, MATLAB® (MathWorks), and LabVIEW® (National Instruments)
- ▶ Integration with EPICS software
- ▶ Rotatable flange for image alignment




Key Camera Features


SOPHIA-XO 2048B - 132



SOPHIA-XO 2048B - 152




SOPHIA-XO 4096B - 154




Large-format, back-illuminated 2048 x 2048 and 4096 x 4096 resolution CCD sensors with >95% peak quantum efficiency in the 5 eV to 30 keV range




Proprietary 4-port readout for low noise and high frame rates



Ultra-high-vacuum, all-metal seal design for deep cooling (ArcTec™) down to -90°C




Flexibility to use air, air+liquid, or liquid cooling




The latest UHV technology with industry-standard CF flange vacuum interface (6 inch or 8 inch)

The most comprehensive software support



Acclaimed LightField for Microsoft Windows 10 (64 bit)



Compatible with EPICS and Python software





PICam software development kit (SDK) for Linux and Microsoft Windows




Exceptional Reliability

Princeton Instruments has been designing low-noise UV / x-ray detectors for more than three decades:

- ▶ Hundreds of cameras being used at leading laboratories around the world
- ▶ Years of trouble-free operation due to uncompromising engineering design and production
- ▶ Complete software ecosystem simplifies image acquisition and processing
- ▶ Continuous innovation to meet evolving requirements and applications



Princeton Instruments SOPHIA-XO soft x-ray cameras are the most advanced CCD detectors yet!




SOPHIA-XO Specifications

Feature	SOPHIA-XO 2048B - 132	SOPHIA-XO 2048B - 152	SOPHIA-XO 4096B - 154
CCD image sensor	e2v CCD42-40; UV-enhanced grade 1; NIMO; back illuminated	e2v CCD230-42; scientific grade 1; AIMO; back illuminated; no AR coating	e2v CCD230-84; scientific grade 1; AIMO; back illuminated; no AR coating
CCD format	2048 x 2048 imaging pixels; 13.5 x 13.5 μm pixels; 100% fill factor	2048 x 2048 imaging pixels; 15.0 x 15.0 μm pixels; 100% fill factor	4096 x 4096 imaging pixels; 15.0 x 15.0 μm pixels; 100% fill factor
Imaging area	27.6 x 27.6 mm	30.7 x 30.7 mm	61.4 x 61.4 mm
Deepest cooling temperature (@ +20°C)	< -90°C (typical) with liquid chiller; < -90°C (typical) with air	< -90°C (typical) with liquid chiller; < -90°C (typical) with air	< -80°C (typical) with liquid chiller; < -60°C (typical) with air
Thermostating precision	$\pm 0.05^\circ\text{C}$		
Dark current (e ⁻ /pixel/sec)	0.0001	0.00025	0.005
Cooling method	Thermoelectric air or liquid cooling		
Full well	Single pixel: 100 ke- (typical)	Single pixel: 150 ke- (typical)	Single pixel: 150 ke- (typical)
ADC speed	8 MHz (4 MHz x 2 ports) 2 MHz (1 MHz x 2 ports) 200 kHz (100 kHz x 2 ports)	16 MHz (4 MHz x 4 ports) 4 MHz (1 MHz x 4 ports) 400 kHz (100 kHz x 4 ports)	
ADC bits	16 bits	16 bits	16 bits 18 bits @ 100 kHz and 1 MHz
System read noise	3.5 e- rms @ 200 kHz 7.0 e- rms @ 2 MHz 19 e- rms @ 8 MHz	3.6 e- rms @ 400 kHz 8.5 e- rms @ 4 MHz 22 e- rms @ 16 MHz	
Readout modes	2-port or 1-port readout; Kinetics; External Sync	4-port, 2-port, or 1-port readout; Kinetics; External Sync	
Nonlinearity	<2% @ 100 kHz		
Software-selectable gains	1, 2, 4 e-/ADU		
Data interface	USB 3.0 (5 m interface cable provided); Optional fiberoptic interface available for remote operation		
I/O signals	Two MCX connectors for programmable frame readout, shutter, trigger in		
Software (optional)	LightField for Microsoft Windows 10 (64 bit); PICam SDK for Microsoft Windows and Linux; EPICS support via automation		
Bake-out temperature	70°C (maximum)		
Vacuum compatibility	10 ⁻⁸ Torr		
Certification	CE		
Operating environment	+5°C to +30°C non-condensing		
Camera head dimensions (L x W x H)	DN100 or 6" industry-standard CF flange: 251.6 mm (9.91") x 129 mm (5.08") x 142.8 mm (5.62")	DN100 or 6" industry-standard CF flange: 251.6 mm (9.91") x 129 mm (5.08") x 142.8 mm (5.62")	DN160 or 8" industry-standard CF flange: 251.6 mm (9.91") x 129 mm (5.08") x 142.8 mm (5.62")

Specifications are subject to change.

Quantum Efficiency Curves



Frame Rates

SOPHIA-XO 2048B - 132

Binning	8 MHz	2 MHz	200 kHz
1 x 1	1.35	0.43	0.05
2 x 2	2.82	1.34	0.18
4 x 4	3.76	2.82	0.65
8 x 8	4.30	3.76	1.82


SOPHIA-XO 2048B - 152

Binning	16 MHz	4 MHz	400 kHz
1 x 1	3.2	0.9	0.09
2 x 2	7.4	2.9	0.33
4 x 4	14.3	7.7	1.05
8 x 8	22.2	15.4	2.9


SOPHIA-XO 4096B - 154

Binning	16 MHz	4 MHz	400 kHz
1 x 1	0.84	0.23	0.024
2 x 2	1.93	0.813	0.082
4 x 4	3.68	2.17	0.258
8 x 8	5.68	4.37	0.700

SOPHIA-XO 4096B



SOPHIA-XO 2048B




LightField[®] Software

The Future of Scientific Imaging and Spectroscopy Software


The combination of LightField and SOPHIA-XO provides researchers with the most advanced and reliable toolset for experimental setup, data acquisition, and post processing:

- ▶ Powerful 64-bit software package includes Microsoft Windows 10 support
- ▶ Complete control of Princeton Instruments cameras and spectrometers
- ▶ Dependable data integrity via automatic saving to disk, time stamping, and retention of both raw and corrected data
- ▶ Full experimental details and system settings are archived and can be reloaded for future experiments ensuring maximum reproducibility
- ▶ For light-sensitive experiments, the user interface offers “low light” and “no light” modes during data acquisition
- ▶ LightField works seamlessly in multi-user facilities, remembering each user’s hardware and software configurations
- ▶ Simple math functions and complex transforms can be applied to live or stored data, includes an easy-to-use editor to create your own formulas
- ▶ Integrated LabVIEW (National Instruments), Python, and MATLAB (MathWorks) support
- ▶ Exports to your favorite file formats, including TIFF, FITS, ASCII, AVI, IGOR, and Origin
- ▶ Demo camera mode allows the user to view all of the settings and parameters associated with any camera without physically connecting the camera
- ▶ Live data processing operations provide real-time evaluation of incoming data to optimize experimental parameters



What is in the box?

SOPHIA-XO cameras can be provided in custom configurations to suit your experiment. Please contact your local Princeton Instruments representative. The most common configurations are listed below:



Optional accessories:

- LightField software
- PICam SDK/API for Linux and Microsoft Windows (provided for free)
- Liquid chiller
- Fiberoptic data extension cable for remote operation from up to 30 m

Other x-ray cameras from Princeton Instruments:

- **PI-MTE & PI-MTE3** – Renowned “in-vacuum” cameras... now in 2k x 2k and 4k x 4k formats
- **PIXIS-XO** – Small-format 1k x 1k and 2k x 2k soft x-ray direct-detection cameras
- **PIXIS-XF & PIXIS-XB** – Fiber-coupled designs with 1k x 1k and 2k x 2k resolutions

Be sure to check out KURO™ (the world’s first back-illuminated sCMOS cameras) and SOPHIA (high-sensitivity CCD cameras) for UV-NIR imaging.

SOPHIA-XO 2048B SOPHIA-XO 4096B

LARGE-FORMAT, DIRECT-DETECTION CCD CAMERAS FOR SOFT X-RAYS



X-ray Image Credits

Front cover, page 4 (upper right), and back cover:
Dr. Jan Luning, SLAC

Page 2:
Yu et al., "Coherent X-ray scattering beamline at port 9C of Pohang Light Source II,"
J. Synchrotron Rad. 21, 264-267 (2014). doi: 10.1107/S1600577513025629

Page 4:
(center) Prof. Jens Biegert and Stephan Teichmann, The Institute of Photonic
Science, Attoscience and Ultrafast Optics, Barcelona, Spain
(lower right) Dr. Jianwei Miao, UCLA

revA4

Copyright © 2018 Princeton Instruments, Inc. All rights reserved. ArcTec and KURO are trademarks and SOPHIA and LightField are registered trademarks of Princeton Instruments, Inc. Microsoft and Windows are registered trademarks of Microsoft Corporation in the United States and other countries. Linux is the registered trademark of Linus Torvalds in the U.S. and other countries. Python is a registered trademark of the Python Software Foundation. All other brand and product names are the trademarks or registered trademarks of their respective owners and manufacturers.

Contact your local Princeton Instruments representative for additional information.

Princeton Instruments - USA

Tel: +1 609.587.9797
pi.info@teledyne.com

China

Tel: +86 157 21 5353 43
pi.info.china@teledyne.com

France

Tel: +33.1.70.38.19.00
EVR@teledyne.com

Germany

Tel: + 49 (0) 89-660 779 3
pi.germany@teledyne.com

Japan

Tel: +81.3.5639.2741
PI.Nippon@teledyne.com

United Kingdom

Tel: +44 (0) 7810 835 719
pi.info@teledyne.com

 **TELEDYNE
PRINCETON INSTRUMENTS**
Everywhereyoulook™