

Single channel detectors are used with IsoPlane, SpectraPro and TriVista spectrometers, enabling efficient operation from UV to IR. These include three photomultiplier tubes, two silicon detectors, plus InGaAs, PbS, InSb, HgCdTe and infrared detectors.

### **Detector Types**

#### **Photomultiplier Tubes**



Photomultiplier tubes are highly sensitive detectors for low-level light in the UV-VIS-NIR range. The tubes are vacuum phototubes which are able to greatly multiply the current produced by any incident light. Photomultiplier tubes are used in applications such as low-light spectroscopy, Raman spectroscopy, astronomy, and particle physics.

#### **Silicon Detectors**



Silicon light detectors are able to convert light into electrical current over the wavelength range 190 – 1100 nm.

They have a long lifetime with high quantum efficiency for a wide range of applications, such as micro-computed tomography and industrial sensors.

# Gallium(Ga), arsenic(As), and indium (In) compounds



Ga, As, and In can be used as part of the detector to cover select wavelengths of the 190 – 2000 nm spectral range. For example, InGaAs detectors are specialised for sensing in the 1000 – 1500 nm range.

## **Family Models**

Model	Description	Available for
ID-441	InGaAs Detector includes pre-amplifier, uncooled	IsoPlane,
	version	SpectraPro, TriVista
ID-441-C	InGaAs Detector includes pre-amplifier, cooling	IsoPlane,
	ready version; Requires 442-1AThermoelectric	SpectraPro, TriVista
ID-442-P	PbS Detector includes pre-amplifier, cooling ready	IsoPlane,
	version. Spectral Range 1.1 – 2.9 µm Requires 442-	SpectraPro, TriVista
	1A thermoelectic cooling for operation. Chopper and	
	lock-in amplifier for operation not included.	
ID-443	InSb Detector, 4mm diameter active area, for use	lsoPlane,
	from 1500-5000nm, LN2 cooled, with housing and	SpectraPro, TriVista
	dewar. Requires chopper and lock-in amplifier for	
<b>DD</b> 400	operation.	
PD-438	Universal PMT detector housing. Mounts to	IsoPlane,
	SpectraPro slit assembly. Requires 1 1/8 side window	SpectraPro, TriVista
PD-439	PMT and high voltage power supply for operation. Universal PMT detector housing. Mounts to	IsoPlane,
FD-439	SpectraPro slit assembly. Requires 1 1/8 side window	SpectraPro, TriVista
	PMT and high voltage power supply for operation.	Specifar 10, minista
PD-471	Integrated PMT detector housing with built-in high	IsoPlane,
	voltage power supply. Mounts to SpectraPro slit	SpectraPro, TriVista
	assembly. Requires 1 1/8 side window PMT and	opooliai io, iii lota
	SpectraHUB for operation.	
PD-473-1	Integrated photon counting detector system for 185-	IsoPlane,
	850 nm. Includes PMT, amplifier/discriminator and	SpectraPro, TriVista
	built-in high voltage power supply. Mounts on	•
	SpectraPro slit assembly. Requires SpectraHUB for	
	operation.	
P1	Photomultiplier Tube, 1 1/8 side window tube,	IsoPlane,
	requires housing and power supply for operation	SpectraPro, TriVista
P2	Photomultiplier Tube, 1 1/8 side window tube,	IsoPlane,
	requires housing and power supply for operation	SpectraPro, TriVista
	HAMAMATSU PN# R928	
P3	Photomultiplier Tube, 1 1/8 side window tube,	IsoPlane,
	requires housing and power supply for operation.	SpectraPro, TriVista
	Range 300 – 1100 nm	
PHV-400	Programmable PMT High Voltage power supply, 0-	IsoPlane,
	1000V	SpectraPro, TriVista
	For use with PMT's and housings	